

PCR Sequencing Kit (SQK-PSK004)

Version: PSK_9072_v1_revA_23May2018
 Last update: 23/05/2018



Flow Cell Number:

DNA Samples:

Before start checklist		
Materials	Consumables	Equipment
<input type="checkbox"/> PCR Sequencing Kit (SQK-PSK004)	<input type="checkbox"/> Agencourt AMPure XP beads	<input type="checkbox"/> Hula mixer (gentle rotator mixer)
<input type="checkbox"/> Flow Cell Priming Kit (EXP-FLP001)	<input type="checkbox"/> NEBNext End repair / dA-tailing Module (E7546)	<input type="checkbox"/> Magnetic rack
	<input type="checkbox"/> NEB Blunt/TA Ligase Master Mix (M0367)	<input type="checkbox"/> Microfuge
	<input type="checkbox"/> Covaris g-TUBE	<input type="checkbox"/> Vortex mixer
	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes	<input type="checkbox"/> Ice bucket with ice
	<input type="checkbox"/> 0.2 ml thin-walled PCR tubes	<input type="checkbox"/> Timer
	<input type="checkbox"/> Nuclease-free water (e.g. ThermoFisher, cat # AM9937)	<input type="checkbox"/> Thermal cycler
	<input type="checkbox"/> Freshly prepared 70% ethanol in nuclease-free water	<input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P200, P1000
	<input type="checkbox"/> LongAmp Taq 2X Master Mix (e.g. NEB M0287)	
	<input type="checkbox"/> 10 mM Tris-HCl pH 8.0 with 50 mM NaCl	
	<input type="checkbox"/> (optional) Exonuclease I (NEB, M0293)	

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Preparing input DNA</p> <p>Prepare the DNA in Nuclease-free water.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Transfer 100 ng DNA into a DNA LoBind tube <input type="checkbox"/> Adjust the volume to 50 µl with Nuclease-free water <input type="checkbox"/> Mix thoroughly by inversion avoiding unwanted shearing <input type="checkbox"/> Spin down briefly in a microfuge <p><input type="checkbox"/> Record the quality, quantity and size of the DNA.</p>	
<p>IMPORTANT</p> <p>Criteria for input DNA</p> <ul style="list-style-type: none"> <input type="checkbox"/> Purity as measured using Nanodrop - OD 260/280 of 1.8 and OD 260/230 of 2.0-2.2 <input type="checkbox"/> Input mass, as measured by Qubit - 100 ng <input type="checkbox"/> No detergents or surfactants in the buffer 	

PCR Sequencing Kit (SQK-PSK004)

Version: PSK_9072_v1_revA_23May2018
 Last update: 23/05/2018



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Check your flow cell</p> <p><input type="checkbox"/> Set up the MinION, flow cell and host computer</p> <p>Once successfully plugged in, you will see a light and hear the fan.</p> <p>Open the MinkNOW GUI from the desktop icon and establish a local or remote connection.</p> <p><input type="checkbox"/> If running a MinION on the same host computer, plug the MinION into the computer.</p> <p><input type="checkbox"/> If running a MinION on a remote computer, first enter the name or IP address of the remote host under Connect to a remote computer (if running from the Connection page), or Connections (if running from the homepage) and click Connect.</p> <p><input type="checkbox"/> Choose the flow cell type from the selector box. Then mark the flow cell as "Selected":</p> <p>Click "Check flow cells" at the bottom of the screen.</p> <p><input type="checkbox"/> R9.4.1 FLO-MIN106</p> <p><input type="checkbox"/> R9.5.1 FLO-MIN107</p> <p><input type="checkbox"/> Click "Start test".</p> <p><input type="checkbox"/> Check the number of active pores available for the experiment, reported in the System History panel when the check is complete.</p>	
<p>Flow cell check complete.</p>	
<p>DNA fragmentation</p> <p>OPTIONAL</p> <p><input type="checkbox"/> Transfer 100 ng genomic DNA in 50 µl to the Covaris g-TUBE.</p> <p>Spin the g-TUBE for 1 minute at RT (Eppendorf 5424; 6000 rpm for 8 kb fragments).</p> <p><input type="checkbox"/> Spin the g-TUBE for 1 minute</p> <p><input type="checkbox"/> Remove and check all the DNA has passed through the g-TUBE</p> <p><input type="checkbox"/> If DNA remains in the upper chamber, spin again for 1 minute at the same speed</p> <p>Invert the g-TUBE and spin again for 1 minute to collect the fragmented DNA.</p> <p><input type="checkbox"/> Remove g-TUBE, invert the tube and replace into the centrifuge</p> <p><input type="checkbox"/> Spin the g-TUBE for 1 minute</p> <p><input type="checkbox"/> Remove and check the DNA has passed into the lower chamber</p> <p><input type="checkbox"/> If DNA remains in the upper chamber, spin again for 1 minute</p> <p><input type="checkbox"/> Remove g-TUBE</p> <p><input type="checkbox"/> Transfer the 50 µl fragmented DNA to a clean 1.5 ml Eppendorf DNA LoBind tube.</p>	
<p>100 ng fragmented DNA in 50 µl is taken into the next step.</p>	
<p>End-prep</p>	

PCR Sequencing Kit (SQK-PSK004)

Version: PSK_9072_v1_revA_23May2018
 Last update: 23/05/2018



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Perform end repair and dA-tailing of fragmented DNA as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 50 µl 100 ng fragmented DNA <input type="checkbox"/> 7 µl Ultra II End-prep reaction buffer <input type="checkbox"/> 3 µl Ultra II End-prep enzyme mix <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Transfer the sample to a 0.2 ml PCR tube, and incubate for 5 minutes at 20 °C and 5 minutes at 65 °C using the thermal cycler.</p> <p><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</p> <p><input type="checkbox"/> Transfer the sample to a 1.5 ml DNA LoBind Eppendorf tube.</p> <p><input type="checkbox"/> Add 60 µl of resuspended AMPure XP beads to the end-prep reaction and mix by pipetting.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual 70% ethanol. Briefly allow to dry.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 16 µl Nuclease-free water. Incubate for 2 minutes at RT.</p> <p><input type="checkbox"/> Pellet beads on magnet until the eluate is clear and colourless.</p> <p><input type="checkbox"/> Remove and retain 16 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <p><input type="checkbox"/> Quantify 1 µl of end-prepped DNA using a Qubit fluorometer.</p> <p><input type="checkbox"/> Repeat the previous step.</p>	
<p>PCR adapter ligation and amplification</p>	
<p>Thaw and prepare the kit reagents as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> PCR Adapter (PCA) at RT <input type="checkbox"/> Whole Genome Primers (WGP) at RT <input type="checkbox"/> Blunt/TA ligation Master Mix on ice <input type="checkbox"/> LongAmp Taq 2x Master Mix on ice 	

PCR Sequencing Kit (SQK-PSK004)

Version: PSK_9072_v1_revA_23May2018
 Last update: 23/05/2018



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Check the contents of each tube are clear of any precipitate and are thoroughly mixed before setting up the reaction.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Mix the contents of each tube by flicking <input type="checkbox"/> Check that there is no precipitate present (DTT in the Blunt/TA Master Mix can sometimes form a precipitate) <input type="checkbox"/> Spin down briefly before accurately pipetting the contents into the reaction <p>Add the reagents in the order given below, mixing by flicking the tube between each sequential addition:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 15 µl End-prepped DNA <input type="checkbox"/> 10 µl PCR Adapters (PCA) <input type="checkbox"/> 25 µl Blunt/TA Ligase Master Mix <ul style="list-style-type: none"> <input type="checkbox"/> Mix gently by flicking the tube, and spin down. <input type="checkbox"/> Incubate the reaction for 10 minutes at RT. <input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing. <input type="checkbox"/> Add 20 µl of resuspended AMPure XP beads to the reaction and mix by pipetting. <input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT. <input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water. <input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant. <input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard. <input type="checkbox"/> Repeat the previous step. <input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. <input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 21 µl Nuclease-free water. Incubate for 2 minutes at RT. <input type="checkbox"/> Pellet beads on magnet until the eluate is clear and colourless. <p>Remove and retain 21 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Remove and retain the eluate which contains the DNA library in a clean 1.5 ml Eppendorf DNA LoBind tube <input type="checkbox"/> Dispose of the pelleted beads <ul style="list-style-type: none"> <input type="checkbox"/> Quantify 1 µl of adapted DNA using a Qubit fluorometer. <input type="checkbox"/> Calculate how much DNA to take forward into the PCR step for a final DNA concentration of 0.2 ng/µl in a 50 µl reaction. <p>Set up the adapted DNA PCR as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Adapter ligated DNA, diluted x µl 0.2 ng/µl <input type="checkbox"/> Nuclease-free water 24-x µl <input type="checkbox"/> Whole Genome Primers (WGP) 1 µl <input type="checkbox"/> LongAmp Taq 2x Master Mix 25 µl 	

PCR Sequencing Kit (SQK-PSK004)

Version: PSK_9072_v1_revA_23May2018
 Last update: 23/05/2018



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p>Amplify using the following cycling conditions:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Initial denaturation 3 mins @ 95 °C (1 cycle) <input type="checkbox"/> Denaturation 15 secs @ 95 °C (14 (b) cycles) <input type="checkbox"/> Annealing 15 secs (a) @ 56 °C (a) (14 (b) cycles) <input type="checkbox"/> Extension 6 mins @ 65 °C (c) (14 (b) cycles) <input type="checkbox"/> Final extension 6 mins @ 65 °C (1 cycle) <input type="checkbox"/> Hold @ 4 °C <p><input type="checkbox"/> Prepare the AMPure XP beads for use; resuspend by vortexing.</p> <p><input type="checkbox"/> Add 30 µl of resuspended AMPure XP beads to the reaction and mix by pipetting.</p> <p><input type="checkbox"/> Incubate on a Hula mixer (rotator mixer) for 5 minutes at RT.</p> <p><input type="checkbox"/> Prepare 500 µl of fresh 70% ethanol in Nuclease-free water.</p> <p><input type="checkbox"/> Spin down the sample and pellet on a magnet. Keep the tube on the magnet, and pipette off the supernatant.</p> <p><input type="checkbox"/> Keep on magnet, wash beads with 200 µl of freshly prepared 70% ethanol without disturbing the pellet. Remove the 70% ethanol using a pipette and discard.</p> <p><input type="checkbox"/> Repeat the previous step.</p> <p><input type="checkbox"/> Spin down and place the tube back on the magnet. Pipette off any residual ethanol. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking.</p> <p><input type="checkbox"/> Remove the tube from the magnetic rack and resuspend pellet in 10 µl of 10 mM Tris.HCl pH 8.0 with 50 mM NaCl. Incubate for 2 minutes at RT.</p> <p><input type="checkbox"/> Pellet beads on magnet until the eluate is clear and colourless.</p> <p>Remove and retain 10 µl of eluate into a clean 1.5 ml Eppendorf DNA LoBind tube.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Remove and retain the eluate which contains the DNA in a clean 1.5 ml Eppendorf DNA LoBind tube <input type="checkbox"/> Dispose of the pelleted beads <p><input type="checkbox"/> Quantify 1 µl of adapted DNA using a Qubit fluorometer.</p> <p><input type="checkbox"/> Make up 50-100 fmol of PCR product to 10 µl in 10 mM Tris-HCl pH 8.0 with 50 mM NaCl.</p>	
<p>Rapid Adapter ligation</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Spin down the Rapid Adapter (RAP) tube, and place it on ice. <input type="checkbox"/> Add 1 µl RAP to the 10 µl amplified DNA library. <input type="checkbox"/> Mix gently by flicking the tube, and spin down. <input type="checkbox"/> Incubate the reaction for 5 minutes at RT. 	

PCR Sequencing Kit (SQK-PSK004)

Version: PSK_9072_v1_revA_23May2018
 Last update: 23/05/2018



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>The prepared library is used for loading into the MinION flow cell. Store the library on ice until ready to load.</p>	
<p>Priming and loading the SpotON flow cell</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Thaw the Sequencing Buffer (SQB), Loading Beads (LB), Flush Tether (FLT) and one tube of Flush Buffer (FLB) at RT before placing the tubes on ice as soon as thawing is complete. <input type="checkbox"/> Mix the Sequencing Buffer (SQB) and Flush Buffer (FLB) tubes by vortexing, spin down and return to ice. <input type="checkbox"/> Spin down the Flush Tether (FLT) tube, mix by pipetting, and return to ice. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Please note that the Sequencing Tether (SQT) tube will NOT be used in this protocol. 	
<ul style="list-style-type: none"> <input type="checkbox"/> Open the lid of the nanopore sequencing device and slide the flow cell's priming port cover clockwise so that the priming port is visible. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Care must be taken when drawing back buffer from the flow cell. The array of pores must be covered by buffer at all times. Removing more than 20-30 µl risks damaging the pores in the array. 	
<p>After opening the priming port, check for small bubble under the cover. Draw back a small volume to remove any bubble (a few µls):</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the priming port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer entering the pipette tip <ul style="list-style-type: none"> <input type="checkbox"/> Prepare the flow cell priming mix: add 30 µl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed Flush Buffer (FLB), and mix by pipetting up and down. <input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for 5 minutes. <input type="checkbox"/> Thoroughly mix the contents of the SQB and LB tubes by pipetting. <p>In a new tube, prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 34 µl Sequencing Buffer (SQB) <input type="checkbox"/> 25.5 µl Loading Beads (LB), mixed immediately before use <input type="checkbox"/> 4.5 µl Nuclease-free water <input type="checkbox"/> 11 µl DNA library 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> The Loading Beads (LB) tube contains a suspension of beads. These beads settle very quickly. It is vital that they are mixed immediately before use. 	
<p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible. <input type="checkbox"/> Load 200 µl of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles. 	

PCR Sequencing Kit (SQK-PSK004)

Version: PSK_9072_v1_revA_23May2018
 Last update: 23/05/2018



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<ul style="list-style-type: none"> <input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading. <input type="checkbox"/> Add 75 µl of sample to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next. <input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION lid. 	
<p>Starting a sequencing run</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Double-click the MinKNOW icon located on the desktop to open the MinKNOW GUI. <input type="checkbox"/> If your MinION was disconnected from the computer, plug it back in. <input type="checkbox"/> Choose the flow cell type from the selector box. Then mark the flow cell as "Selected". <input type="checkbox"/> Click the "New Experiment" button at the bottom left of the GUI. <p>Select the appropriate protocol script</p> <ul style="list-style-type: none"> <input type="checkbox"/> Experiment type: Choose Sequencing Run under "Choose Operation" <input type="checkbox"/> Flow Cell product code: Choose the Flow Cell type under "Flow cell product code" <input type="checkbox"/> Sequencing kit: Choose SQK-PSK004 under Sequencing Kit <input type="checkbox"/> Choose whether or not live basecalling is enabled <input type="checkbox"/> The most appropriate script will appear in the drop-down menu. <p><input type="checkbox"/> Click "Begin Experiment".</p> <p>Allow the script to run to completion.</p> <ul style="list-style-type: none"> <input type="checkbox"/> The MinKNOW Experiment page will indicate the progression of the script; this can be accessed through the "Experiment" tab that will appear at the top right of the screen <input type="checkbox"/> Monitor messages in the Message panel in the MinKNOW GUI <p>The basecalled read files are stored in : \data\reads</p>	
<p>Progression of MinKNOW protocol script</p>	
<p>The running experiment screen</p> <p>Experiment summary information</p> <p>Check the number of active pores reported in the MUX scan are similar (within 10-15%) to those reported at the end of the Platform QC</p> <ul style="list-style-type: none"> <input type="checkbox"/> If there is a significant reduction in the numbers, restart MinKNOW. <input type="checkbox"/> If the numbers are still significantly different, close down the host computer and reboot. <input type="checkbox"/> When the numbers are similar to those reported at the end of the Platform QC, restart the experiment on the Connection page. There is no need to load any additional library after restart. <input type="checkbox"/> Stopping the experiment is achieved by clicking "Stop experiment" button at the top of the screen. <p><input type="checkbox"/> Check the temperature is approximately 34° C.</p>	

Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Check pore occupancy in the channel panel at the top of the experimental view.</p> <ul style="list-style-type: none"> <input type="checkbox"/> A good library will be indicated by a higher proportion of light green channels in Sequencing than are in Pore. The combination of Sequencing and Pore indicates the number of active pores at any point in time. A low proportion of Sequencing channels will reduce the throughput of the run. <input type="checkbox"/> Recovering indicates channels that may become available for sequencing again. A high proportion of this may indicate additional clean up steps are required during your library preparation. <input type="checkbox"/> Inactive indicates channels that are no longer available for sequencing. A high proportion of these as soon as the run begins may indicate an osmotic imbalance. <input type="checkbox"/> Unclassified are channels that have not yet been assigned one of the above classifications <p><input type="checkbox"/> Monitor the pore occupancy</p> <p>Duty time plots</p> <p><input type="checkbox"/> Monitor the development of the read length histogram.</p> <p>Trace viewer</p>	
<p>Further analysis with EPI2ME (optional)</p>	
<p>OPTIONAL</p> <ul style="list-style-type: none"> <input type="checkbox"/> Open the Desktop Agent using the desktop shortcut. <input type="checkbox"/> Click on the New Workflow tab in the Desktop Agent and select the workflow to be used in the analysis. <p>Select the workflow parameters.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Select the quality score cut-off (this defaults to 7 unless changed) <input type="checkbox"/> Select whether the reads should be demultiplexed by barcode <input type="checkbox"/> If you are working with human data, please tick "Yes" in answer to "Is the data you are about to upload a whole or partial human genome?", and confirm that you have consent from the subject to upload the data. <ul style="list-style-type: none"> <input type="checkbox"/> Check the correct settings are selected in the Desktop Agent. <input type="checkbox"/> Click "Start Run" to start data analysis. <input type="checkbox"/> Follow the progression of upload and download of read files in the Desktop Agent. <p>Click on VIEW REPORT.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Click on VIEW REPORT to navigate to the EPI2ME website, this can be done at any point during data exchange <input type="checkbox"/> Return to the Desktop Agent to see progression of the exchange <ul style="list-style-type: none"> <input type="checkbox"/> When the upload and download numbers are the same, the data exchange is complete. The processed reads will be in downloads folder in the selected location on the host computer. 	
<p>Close down MinKNOW and the Desktop Agent</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Quit Desktop Agent using the close x. <input type="checkbox"/> Quit MinKNOW by closing down the web GUI. 	

PCR Sequencing Kit (SQK-PSK004)

Version: PSK_9072_v1_revA_23May2018
Last update: 23/05/2018



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<input type="checkbox"/> Disconnect the MinION.	
Prepare the flow cell for re-use or return to Oxford Nanopore.	
<input type="checkbox"/> If you would like to reuse the flow cell, follow the Wash Kit instructions and store the washed flow cell at 2-8 °C, OR <input type="checkbox"/> Follow the returns procedure by washing out the MinION Flow Cell ready to send back to Oxford Nanopore.	