

Rapid Lambda Control Experiment (SQK-RAD004)

Version: RSE_9045_v1_revA_17Nov2017
 Last update: 17/11/2017



Flow Cell Number:

DNA Samples:

Before start checklist		
<input type="checkbox"/> Rapid Sequencing Kit (SQK-RAD004)	<input type="checkbox"/> Timer	<input type="checkbox"/> 1.5 ml Eppendorf DNA LoBind tubes
<input type="checkbox"/> Flow Cell Priming Kit (EXP-FLP001)	<input type="checkbox"/> Method of heating to 80° C for 1 minute	<input type="checkbox"/> 0.2 ml thin-walled PCR tubes
<input type="checkbox"/> Microfuge	<input type="checkbox"/> Pipettes and pipette tips P2, P10, P20, P100, P1000	<input type="checkbox"/> Nuclease-free water

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Check your flow cell</p> <p><input type="checkbox"/> Set up the MinION, Flow Cell and host computer</p> <p>Once successfully plugged in, you will see a light and hear the fan.</p> <p>Open the MinkNOW GUI from the desktop icon and establish a local or remote connection.</p> <ul style="list-style-type: none"> <input type="checkbox"/> If running a MinION on the same host computer, plug the MinION into the computer. When the connection name appears under the Local tab, click Connect. <input type="checkbox"/> If running a MinION on a remote computer, first enter the name or IP address of the remote host under Remote and click Connect. <input type="checkbox"/> Plug a MinION and Flow Cell into the remote computer; the connection IDs will be displayed under MinION Connection and Flowcell Connection. <p>Enter the SampleID and FlowcellID being used, and click Submit.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Once a MinION and Flow Cell are connected, a Label Experiment dialogue box appears. <input type="checkbox"/> Click into the Sample ID box and name your sample using free text in alphanumeric format only, deleting any default Sample_ID that is present. Warning: SampleID should not contain any personally identifiable information. <input type="checkbox"/> Click into the FlowcellID box and enter the Flow Cell ID, which is the code found on a sticker on the top side of a Flow Cell. <p><input type="checkbox"/> Select the Platform QC script under Choose Operation, and start the script using the Execute button.</p> <p><input type="checkbox"/> Check the number of active pores available for the experiment, reported in the message panel or in notifications when the check is complete.</p>	
Flow cell check complete.	
Library preparation	
DNA tagmentation	

Rapid Lambda Control Experiment (SQK-RAD004)

Version: RSE_9045_v1_revA_17Nov2017
 Last update: 17/11/2017



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p>Thaw kit components at RT, spin down briefly using a microfuge and mix by pipetting as indicated by the table below:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Lambda DNA (LMD): thaw at RT, briefly spin down, mix well by pipetting <input type="checkbox"/> Fragmentation Mix (FRA): not frozen, briefly spin down, mix well by pipetting <input type="checkbox"/> Rapid Adapter (RAP): not frozen, briefly spin down, mix well by pipetting <input type="checkbox"/> Sequencing Buffer (SQB): thaw at RT, briefly spin down, mix well by pipetting* <input type="checkbox"/> Loading Beads (LB): thaw at RT, briefly spin down, mix by pipetting or vortexing immediately before use <input type="checkbox"/> Flush Buffer (FLB): thaw at RT, briefly spin down, mix well by pipetting* <input type="checkbox"/> Flush Tether (FLT): thaw at RT, briefly spin down, mix well by pipetting <p><input type="checkbox"/> Once thawed, keep all the kit components on ice.</p> <p>In a 0.2 ml thin-walled PCR tube, mix the following:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 7.5 µl Lambda DNA <input type="checkbox"/> 2.5 µl FRA <p><input type="checkbox"/> Mix gently by flicking the tube, and spin down.</p> <p><input type="checkbox"/> Incubate the tube at 30° C for 1 minute and then at 80° C for 1 minute.</p>	
<p>375 ng tagmented Lambda DNA in 10 µl is taken into the next step.</p>	
<p>Adapter attachment</p> <ul style="list-style-type: none"> <input type="checkbox"/> Add 1 µl of RAP to the tube. <input type="checkbox"/> Mix gently by flicking the tube, and spin down. <input type="checkbox"/> Incubate the reaction for 5 minutes at RT. 	
<p>The prepared DNA library is used for loading into the flow cell. Store the library on ice until ready to load.</p>	
<p>Priming and loading the SpotON Flow Cell</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Open the lid of the nanopore sequencing device and slide the flow cell's priming port cover clockwise so that the priming port is visible. 	
<p>IMPORTANT</p> <ul style="list-style-type: none"> <input type="checkbox"/> Care must be taken when drawing back buffer from the flow cell. The array of pores must be covered by buffer at all times. Removing more than 20-30 µl risks damaging the pores in the array. 	
<p>After opening the priming port, check for small bubble under the cover. Draw back a small volume to remove any bubble (a few µl):</p> <ul style="list-style-type: none"> <input type="checkbox"/> Set a P1000 pipette to 200 µl <input type="checkbox"/> Insert the tip into the priming port <input type="checkbox"/> Turn the wheel until the dial shows 220-230 µl, or until you can see a small volume of buffer entering the pipette tip <p><input type="checkbox"/> Prepare the flow cell priming mix: add 30 µl of thawed and mixed Flush Tether (FLT) directly to the tube of thawed and mixed Flush Buffer (FLB), and mix by pipetting up and down.</p>	

Rapid Lambda Control Experiment (SQK-RAD004)

Version: RSE_9045_v1_revA_17Nov2017
 Last update: 17/11/2017



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
<p><input type="checkbox"/> Load 800 µl of the priming mix into the flow cell via the priming port, avoiding the introduction of air bubbles. Wait for 5 minutes.</p> <p><input type="checkbox"/> Thoroughly mix the contents of the SQB and LB tubes by pipetting.</p> <p>In a new tube, prepare the library for loading as follows:</p> <ul style="list-style-type: none"> <input type="checkbox"/> 34 µl Sequencing Buffer (SQB) <input type="checkbox"/> 25.5 µl Loading Beads (LB) <input type="checkbox"/> 4.5 µl Nuclease-free water <input type="checkbox"/> 11 µl DNA library <p>Complete the flow cell priming:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Gently lift the SpotON sample port cover to make the SpotON sample port accessible. <input type="checkbox"/> Load 200 µl of the priming mix into the flow cell via the priming port (not the SpotON sample port), avoiding the introduction of air bubbles. <input type="checkbox"/> Mix the prepared library gently by pipetting up and down just prior to loading. <input type="checkbox"/> Add 75 µl of sample to the flow cell via the SpotON sample port in a dropwise fashion. Ensure each drop flows into the port before adding the next. <input type="checkbox"/> Gently replace the SpotON sample port cover, making sure the bung enters the SpotON port, close the priming port and replace the MinION lid. 	
<p>Onward analysis of MinKNOW basecalled data</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Open the Desktop Agent using the desktop shortcut. <input type="checkbox"/> Click on the New Workflow tab in the Desktop Agent and select the Control Experiment workflow. <input type="checkbox"/> Check the correct settings are selected in the Desktop Agent. <input type="checkbox"/> Click "Start Run" to start data analysis. <input type="checkbox"/> Follow the progression of upload and download of read files in the Desktop Agent, along with network speed. <p>Click on VIEW REPORT.</p> <ul style="list-style-type: none"> <input type="checkbox"/> Click on VIEW REPORT to navigate to the Metrichor website, this can be done at any point during data exchange <input type="checkbox"/> Return to the Desktop Agent to see progression of the exchange 	
<p>Close down MinKNOW and the Desktop Agent</p>	
<ul style="list-style-type: none"> <input type="checkbox"/> Quit Desktop Agent using the close x. <input type="checkbox"/> Quit MinKNOW by closing down the web GUI. <input type="checkbox"/> Disconnect the MinION. 	

Rapid Lambda Control Experiment (SQK-RAD004)

Version: RSE_9045_v1_revA_17Nov2017
Last update: 17/11/2017



Flow Cell Number:

DNA Samples:

INSTRUCTIONS	NOTES/OBSERVATIONS
Prepare the flow cell for re-use or return to Oxford Nanopore.	
<input type="checkbox"/> If you would like to reuse the flow cell, follow the Wash Kit instructions and store the washed flow cell at 2-8 °C, OR <input type="checkbox"/> Follow the returns procedure by washing out the MinION Flow Cell ready to send back to Oxford Nanopore.	